



Transducer Holder

CLINICAL ADVANTAGES

| Features | Benefits | Clinical Advantages |
|--------------------------|---|--|
| Secure Platform | <ul style="list-style-type: none"> • Secures and organizes up to three transducers • Platform is soft and padded • Soft elasticized strap to secure to patient's upper arm | <ul style="list-style-type: none"> • Easy transducer positioning at phlebostatic axis enhances leveling accuracy • Monitoring at phlebostatic axis is maintained during patient repositioning or transport |
| 3 Holding Tabs | <ul style="list-style-type: none"> • Secures 1, 2 or 3 transducers • Color coded • Easy to use Velcro® fasteners | <ul style="list-style-type: none"> • Application is easily identified, colors coordinate CVP, Arterial, and PA pressure lines • Facilitates user learning curve • Less manipulation of transducers and pressure tubing |
| Easy to Apply | <ul style="list-style-type: none"> • Affixed to patient extremity • Can be used over BP Cuff or sleeve | <ul style="list-style-type: none"> • Reduces the time and trouble of re-leveling with patient or bed position changes • Eliminates taping individual transducers to the patient's chest or arm or mounting on bedside pole • Saves nursing time |
| Convenient System | <ul style="list-style-type: none"> • The holder is designed to secure and organize the transducers at mid-chest area | <ul style="list-style-type: none"> • Allows for a complete securement system for radial artery blood pressure monitoring when used in conjunction with the Dale Bendable ArmBoard. |
| Latex-Free | <ul style="list-style-type: none"> • Eliminates patient and caregiver latex reactions | <ul style="list-style-type: none"> • Prevents latex hypersensitivity |

Velcro is a registered trademark of Velcro USA, Inc.

DALE MEDICAL PRODUCTS, INC., 7 CROSS STREET, PO BOX 1556, PLAINVILLE, MA 02762-0056 USA
 TEL: 800-343-3980 INT'L: 508-695-9316 FAX: 508-695-6587 www.dalemed.com